


Comparison of monitoring systems functions for dairy cows



e-stado® works for you whole year round.

Activity meters work for you 2-3 months per year.

*integration with milking systems

e-mail: info@e-stado.net contact us for details: www.e-stado.net



Dairy cow monitoring

Easy

Ear biosensor E2



Functions:

- Reproduction
- Milk*

Opti

Ear biosensor E3



Functions:

- Reproduction
- Nutrition
- Activity
- Milk*

Pro

Ear biosensor E4



Functions:

- Reproduction
- Nutrition
- Activity
- Ear temperature
- Individual summer heat stress
- Milk*

Calve

Tail biosensor T4



Functions:

- Calving
- Postpartum retention milk fever
- Milk*

*integration with milking systems

Options for all biosensors:

- Cowshed temperature and humidity (TH2 sensor needed)
- Cowshed photoperiod (TL2 sensor needed)

contact us for details: e-mail: info@e-stado.net www.e-stado.net